Bittium

Service HSAT de la plateforme MedicalSuite™

Guide clinique

Publié par

Bittium Biosignals Ltd.
Pioneerinkatu 6
70800 Kuopio
Finlande
Tél. : +358 40 344 2000
www.bittium.com

Avis juridique

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The screenshot shows the Bittium Recordings page. At the top, there is a navigation bar with the Bittium logo, 'Dashboard', 'Recordings', and user information 'EN - Sales Technician Technician'. Below the navigation bar, there is a search bar with a magnifying glass icon and the text 'Search', a 'Create recording' button, and two dropdown menus for '10' and 'All'. The main content is a table with the following columns: Patient ID, Status, Last data transfer, Last update, Organization, Created, Type, Device, Assignment, and Action. The table contains 10 rows of data, each representing a recording. The status of each recording is indicated by a colored dot: green for 'Active' and red for 'Closed'. The 'Action' column contains icons for downloading and deleting records.

Patient ID	Status	Last data transfer	Last update	Organization	Created	Type	Device	Assignment	Action
KK02	Active	28.10.2021 15:25	28.10.2021 15:46	Sleep sales center	28.10.2021 15:23	Bittium HSAT	0000568	Sales diag group 1	Download
Pat123456	Closed	19.10.2021 11:41	19.10.2021 12:12	Sleep sales center	19.10.2021 11:40	Bittium HSAT	0000673	Sales diag group 1	Delete Download
PID0000b	Closed	18.10.2021 17:48	18.10.2021 18:16	Sleep sales center	18.10.2021 17:47	Bittium HSAT	0000673	Sales diag group 1	Delete Download
PID0000a	Closed	18.10.2021 16:19	18.10.2021 16:15	Sleep sales center	18.10.2021 16:18	Bittium HSAT	0000536	Sales diag group 1	Delete Download
PID00009	Active	16.10.2021 08:30	09.11.2021 10:16	Sleep sales center	16.10.2021 08:29	Bittium HSAT	0000536	Sales diag group 1	Download
PID00008	Active	16.10.2021 08:25	16.10.2021 09:04	Sleep sales center	16.10.2021 08:25	Bittium HSAT	0000568	Sales diag group 1	Download
Test008	Closed	15.10.2021 19:01	18.10.2021 15:26	Sleep sales center	15.10.2021 19:00	Bittium HSAT	0000673	Sales diag group 1	Delete Download
TEST007	Closed	15.10.2021 18:46	15.10.2021 19:14	Sleep sales center	15.10.2021 18:43	Bittium HSAT	0000673	Sales diag group 1	Delete Download
PID00007	Active	No date	14.10.2021 17:22	Sleep sales center	14.10.2021 17:18	Bittium HSAT	1234567	Sales diag group 1	Download
23_45_08	Closed	01.10.2021 16:17	14.10.2021 17:19	Sleep sales center	01.10.2021 16:16	Bittium HSAT	0123456	Sales diag group 1	Delete Download

Figure 13 Vue Enregistrements de MedicalSuite Center

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Patient ID: Pat123456 Bittium HSAT

Approved

Organization	Sleep sales center	Created	19.10.2021 11:40	Last data transfer	19.10.2021 11:41
Status	● Closed	Started	18.10.2021 23:30	Last update	19.10.2021 12:12
Assignment	Sales diag group 1 ▾	Ready for analysis	No date	Duration	1 / 1 day(s)

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Cancel Save

Patient details

Patient ID	Pat123456	Age (years)	55
	<small>Patient ID must be at least 4 characters long and lower/uppercase letters, numbers and underscores are allowed. If Patient ID is left blank the system will generate it (for example PID0000A).</small>	Height (cm)	175
		Weight (kg)	80
		BMI	26.12
		Gender	Female ▾

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
A screenshot of the Bittium MedicalSuite Center interface. At the top, there is a navigation bar with 'Bittium', 'Dashboard', and 'Recordings' tabs. On the right, it shows 'EN' and a user profile for 'Sales Technician'. Below this, the main heading is 'Create recording: Bittium HSAT'. There are four tabs: 'Patient details*' (selected), 'Anamneses', 'Recording details*', and 'Attachments'. Navigation buttons 'Back' and 'Next' are visible, along with 'Cancel' and 'Save' buttons. The 'Patient details' section contains several input fields: 'Patient ID' with a note that it must be at least 4 characters long and can contain lowercase/uppercase letters, numbers, and underscores; 'Age (years)*', 'Height (cm)*', 'Weight (kg)*', 'BMI*', and 'Gender*'. Each field has a corresponding input line.

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The screenshot displays the Bittium web application interface for creating a recording. The top navigation bar includes the Bittium logo, 'Dashboard', 'Recordings', and a user profile for 'Sales Technician'. The main heading is 'Create recording: Bittium HSAT'. Below this, there are tabs for 'Patient details', 'Anamneses', 'Recording details', and 'Attachments'. Navigation buttons for '<< Back', 'Next >>', 'Cancel', and 'Save' are present. The 'Anamnesis' tab is active, showing two sections: 'Medication' and 'Medical history'. The 'Medication' section has a text input field and a 'Medication during the recording' section with radio buttons for 'No' and 'Yes'. The 'Medical history' section includes a question about family members with 'No' and 'Yes' radio buttons, a text input field, and a list of conditions with 'No' and 'Yes' radio buttons: Hypertension, Diabetes, Asthma, Hypothyroidism, Stuffy nose, Enlarged tonsils, and Depression. Below this is the 'Life' section with 'Smoking' and 'Professional driver' radio buttons, and a 'Driving license class' text input field.

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Chance of dozing off or falling asleep in the following situations

Sitting and reading	<input type="radio"/> Never	<input type="radio"/> Slight	<input type="radio"/> Moderate	<input type="radio"/> High
Watching TV	<input type="radio"/> Never	<input type="radio"/> Slight	<input type="radio"/> Moderate	<input type="radio"/> High
Sitting, inactive in a public place (e.g. a theatre or a meeting)	<input type="radio"/> Never	<input type="radio"/> Slight	<input type="radio"/> Moderate	<input type="radio"/> High
As a passenger in a car for an hour without a break	<input type="radio"/> Never	<input type="radio"/> Slight	<input type="radio"/> Moderate	<input type="radio"/> High
Lying down to rest in the afternoon when circumstances permit	<input type="radio"/> Never	<input type="radio"/> Slight	<input type="radio"/> Moderate	<input type="radio"/> High
Sitting and talking to someone	<input type="radio"/> Never	<input type="radio"/> Slight	<input type="radio"/> Moderate	<input type="radio"/> High
Sitting quietly after a lunch without alcohol	<input type="radio"/> Never	<input type="radio"/> Slight	<input type="radio"/> Moderate	<input type="radio"/> High
In a car, while stopped for a few minutes in the traffic	<input type="radio"/> Never	<input type="radio"/> Slight	<input type="radio"/> Moderate	<input type="radio"/> High

Sleep history

Morning headache	<input type="radio"/> Never	<input type="radio"/> Rarely	<input type="radio"/> Sometimes	<input type="radio"/> Often	<input type="radio"/> Almost always
Dry mouth in the morning	<input type="radio"/> Never	<input type="radio"/> Rarely	<input type="radio"/> Sometimes	<input type="radio"/> Often	<input type="radio"/> Almost always
Daytime fatigue or snoring	<input type="radio"/> Never	<input type="radio"/> Rarely	<input type="radio"/> Sometimes	<input type="radio"/> Often	<input type="radio"/> Almost always
Daytime napping	<input type="radio"/> Never	<input type="radio"/> Rarely	<input type="radio"/> Sometimes	<input type="radio"/> Often	<input type="radio"/> Almost always
Heartburn	<input type="radio"/> Never	<input type="radio"/> Rarely	<input type="radio"/> Sometimes	<input type="radio"/> Often	<input type="radio"/> Almost always
Loss of concentration	<input type="radio"/> Never	<input type="radio"/> Rarely	<input type="radio"/> Sometimes	<input type="radio"/> Often	<input type="radio"/> Almost always
Memory loss	<input type="radio"/> Never	<input type="radio"/> Rarely	<input type="radio"/> Sometimes	<input type="radio"/> Often	<input type="radio"/> Almost always
Cramping at night	<input type="radio"/> Never	<input type="radio"/> Rarely	<input type="radio"/> Sometimes	<input type="radio"/> Often	<input type="radio"/> Almost always
Difficulties falling asleep	<input type="radio"/> Never	<input type="radio"/> Rarely	<input type="radio"/> Sometimes	<input type="radio"/> Often	<input type="radio"/> Almost always
Breathing breaks while sleeping	<input type="radio"/> Never	<input type="radio"/> Rarely	<input type="radio"/> Sometimes	<input type="radio"/> Often	<input type="radio"/> Almost always
Pain during night	<input type="radio"/> Never	<input type="radio"/> Rarely	<input type="radio"/> Sometimes	<input type="radio"/> Often	<input type="radio"/> Almost always
Snoring or snorting during sleep	<input type="radio"/> Never	<input type="radio"/> Rarely	<input type="radio"/> Sometimes	<input type="radio"/> Often	<input type="radio"/> Almost always
Waking up during the night or restless sleep	<input type="radio"/> Never	<input type="radio"/> Rarely	<input type="radio"/> Sometimes	<input type="radio"/> Often	<input type="radio"/> Almost always
A strong urge to move your feet or an uncomfortable sensation in your feet related to this urge during sleep	<input type="radio"/> Never	<input type="radio"/> Rarely	<input type="radio"/> Sometimes	<input type="radio"/> Often	<input type="radio"/> Almost always
Urinary urgency during sleep	<input type="radio"/> Never	<input type="radio"/> Rarely	<input type="radio"/> Sometimes	<input type="radio"/> Often	<input type="radio"/> Almost always
If yes, how many times per night?	<input type="radio"/> Once	<input type="radio"/> Twice	<input type="radio"/> More often		
Average amount of sleep	_____ Hours				

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Create recording: Bittium HSAT

Patient details* Anamneses **Recording details** Attachments

<< Back Next >>

Recording details

Organization ServiceProvider Referring physician

Requested report language English

Reason for recording

Device information

Recorder S/N (21) 19 XXXXXXXX

See the side of the device. **Please add numbers after "19"**.

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Create recording: Bittium HSAT

Patient details Anamneses Recording details* Attachments

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Patient ID: **PID0a74a** Bittium HSAT

Ready for analysis ⋮

Organization	Sleep test center	Created	29.12.2021 12:06	Last data transfer	29.12.2021 12:20
Status	● Active ▼	Started	09.09.2021 23:30	Last update	29.12.2021 12:21
Assignment	None ▼	Ready for analysis	No date	Duration	2 / 2 day(s)

[Recording data](#)
[Patient details](#)
[Anamneses](#)
Patient diary
[Recording details](#)
[Attachments](#)

Cancel
Save

Patient diary - 09.09.2021, 23:30 (06h 59m)

Please, ask the patient to answer the following questions about the recording.

When did you go to bed? <input style="width: 100%;" type="text" value="hh:mm"/>	In your estimation, when did you fall asleep? <input style="width: 100%;" type="text" value="hh:mm"/>
What time did you wake up in the morning? <input style="width: 100%;" type="text" value="hh:mm"/>	In your estimation, how many hours did you sleep? <input style="width: 100%;" type="text" value="hh:mm"/>
How did you sleep? <input style="width: 100%;" type="text"/>	Did you get out of bed during the recording? <input style="width: 100%;" type="text"/>
How did you sleep compared to usual? <input style="width: 100%;" type="text"/>	

Other comments

If the patient woke up during the recording, please write down the events.

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Status	Uploading	Created	01.04.2019 15:34	Last update	01.04.2019 15:34
Assignment	Mabel Hicks	Started	---	Duration	---

Recording data Patient details Anamneses Patient diary Recording details Attachments

Download
 Include all attachments

Cancel Save

Recording data

UPLOADS

Do not close this browser tab

Patient ID: 123455677888

- Uploading
123455677886.edf
123455677887.edf
- Data preparation

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Recording data* Patient details Anamneses Patient diary Recording details Attachments

Download
 Include all attachments

Cancel Save

Recording data

+ Add file or drag here

123455677886.edf ✘
The recording file is not valid, please choose another file.

document_2.edf ✘
The used recorder 1234567 doesn't match the selected recorder 9876543. Please make sure that the file is correct or change the recorder.

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The screenshot displays the Bittium MedicalSuite interface. At the top, a metadata table shows details for 'Clinic A', including creation and update dates. Below this, a navigation bar highlights 'Recording data'. A 'Download' section is visible with a 'Save' button. The main area features a table of recording data with columns for Day no., Start time, Duration, Status, Quality, and Notes. The second row is highlighted, showing a recording from 11.04.2019 with a 'Good' quality status. A red box highlights the 'Preparing' status icon in the table. To the right, a modal window titled 'UPLOADS' shows the upload progress for Patient ID: 12345677888, with one file successfully uploaded and another in the 'Preparing data' stage.

Organization	Clinic A	Created	01.04.2019 15:34	Last data transfer	13.04.2019 12:00
Status	Preparing	Started	10.04.2019 12:00	Last update	13.04.2019 12:00
Assignment	None	Ended	No date	Duration	2 days

Day no.	Start time	Duration	Status	Quality	Notes
1 / 2	10.04.2019 22:00	08h 05m	Open	Fair	Add notes
2 / 2	11.04.2019 22:00	07h 50m	Open	Good	Add note

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EN Sales Technician

Patient ID: PID0a74a Bittium HSAT Ready for analysis

Organization	Sleep test center	Created	29.12.2021 11:22	Last data transfer	29.12.2021 11:32
Status	Active	Started	16.12.2021 02:44	Last update	29.12.2021 11:33
Assignment	None	Ready for analysis	No date	Duration	2 / 2 day(s)

Recording data Patient details Anamneses Patient diary Recording details Attachments

Download Include all attachments Cancel Save

Recording data

+ Add file or drag here

<input type="checkbox"/>	Day No.	Start date of file	Duration	Status	Quality	Notes
<input type="checkbox"/>	1/1	16.12.2021 00:44	08h 10m	Open	Good	Add notes

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The screenshot shows a software interface with a navigation bar at the top containing tabs: Recording data, Patient details, Anamneses, Patient diary, Recording details, and Attachments. Below the navigation bar, there are buttons for 'Download', 'Include all attachments', 'Cancel', and 'Save'. The main content area is titled 'Recording data' and contains a table with columns: Day No., Start date of file, Duration, Status, Quality, and Notes. The table has three rows of data. The 'Status' column for the first row is 'Done' (red dot), the second row is 'Open' (yellow dot), and the third row is 'Unusable' (purple dot). A red box highlights the 'Status' dropdown menu for the second row, which is open and shows options: 'Open', 'Pending', 'Done', 'Unusable', and 'Unusable'. To the right of the table, there is a section titled 'Add notes' with a sub-section 'Marked as "Unusable"' containing a list of reasons: 'Corrupted recording data', 'Not enough analyzable data due to noisy signal or unattached recorder', 'Not enough sleep to write a statement', and 'Other (ts)'. There are also icons for adding notes and deleting rows.

Day No.	Start date of file	Duration	Status	Quality	Notes
1/3	14.12.2021 23:42	04h 30m	Done	Fair	Add notes
2/3	14.12.2021 23:42	05h 00m	Open	Good	Add notes
3/3	14.12.2021 23:42	00h 15m	Unusable	Poor	Add notes

Marked as "Unusable"

- Corrupted recording data
- Not enough analyzable data due to noisy signal or unattached recorder
- Not enough sleep to write a statement
- Other (ts)

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Unusable recording data

Patient ID: [patient_id], 10/04/2019, 22:00 (08h 20m)

Why the recording data of this night is unusable?

Corrupted recording data

Not enough analyzable data due to noisy signal or unattached recorder

Not enough sleep to write the statement

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The screenshot displays the Bittium user interface for a patient's recording. At the top, there is a navigation bar with 'Back', 'Bittium', 'Dashboard', and 'Recordings'. The user is identified as 'Sales Technician'. The patient ID is 'PID0a74a' and the recording is for 'Bittium HSAT'. A red box highlights the 'Ready for analysis' button in the top right corner. Below this, a table shows recording details: Organization (Sleep test clinic A), Created (24.03.2022 15:55), Last data transfer (24.03.2022 15:58), Status (Active), Started (14.12.2021 23:42), Last update (05.05.2022 13:15), Assignment, Ready for analysis (No date), and Duration (2 / 2 day(s)). A dropdown menu for the status is open, showing 'Active' and 'Ended' options, with 'Ended' highlighted by a red box. Below the table, there are tabs for 'Recording data', 'Patient details', 'Anamneses', 'Patient diary', 'Recording details', and 'Attachments'. A 'Download' section includes a checkbox for 'Include all attachments' and 'Cancel'/'Save' buttons. At the bottom, the 'Recording data' section shows a table with columns for Day No., Start date of file, Duration, Status, Quality, and Notes. One recording is listed with Day No. 1/2, Start date 14.12.2021 23:42, Duration 04h 30m, Status Open, and Quality Fair.

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Ready for analysis

Patient ID **PID0a870_013696**

Are you sure this recording is ready for analysis?

- Please make sure the patient diary is filled or it is attached to "Other documents" section
- You will no longer be able to upload new recording data to it.

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Diagnostic report needed


Recording	Assignment	Ready for analysis	Duration	
● PID0a74a	Test diag group 1 ▾	24.11.2021 09:10	1 day(s)	Analyze 

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The screenshot shows the Bittium interface for Patient ID: PID00298. The top navigation bar includes 'Back', 'Bittium', 'Dashboard', 'Recordings', 'EN', a help icon, and the user profile 'Technician 1'. The patient information section shows 'Bittium HSAT' and a 'Ready for analysis' button. Below this is a table with the following data:

Organization	DemoClinic1	Created	30.05.2023 13:56	Last data transfer	30.05.2023 14:01
Status	Active	Started	23.11.2021 23:42	Last update	25.01.2024 14:08
Assignment	HSAT Analysis	Ready for analysis	No date	Duration	1 / 1 day(s)

A dropdown menu is open for the 'Assignment' field, showing a search bar and a list of groups: 'Cardiscope Group', 'Diagnostic group', 'Diagnostic group (autoreports)', 'HRV group (autoreports)', and 'HSAT Analysis' (selected). Below the groups, it shows 'Users in HSAT Analysis' with 'Dr. Hanna HSAT' listed. The background shows a 'Patient diary' tab and 'Recording details' and 'Attachments' sections. There are 'Cancel' and 'Save' buttons at the bottom right of the assignment field.

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Diagnostic report needed

All	Not assigned				
Recording	Assignment	Ready for analysis	Duration		
● PatientZero	Test diag group 1	24.11.2021 09:10	1 day(s)	Analyze	↓

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1 of 3 Automatic Zoom

Patient

Patient ID:	23-30-48	Height:	185cm	Weight:	85kg
BMI:	24.8	Age:	58 years	Gender:	Male

Summary

Apnea-hypopnea index: 8.7/h
(Supine: 15.4/h, Non-supine: 2.1/h)

Oxygen Desaturation Index (ODI3): 8.7/h

AHI

Normal 5 Mild 15 Moderate 30 Severe

ODI3

Recording date	15/10/2021	Recording: TRT	08h 36m (23:30 - 08:07)
Evaluation period	08h 16m (23:43 - 08:00)		

Apnea-hypopnea index (AHI) displayed here is the total number of respiratory events (Apneas + Hypopneas) divided by Evaluation period. Evaluation period was determined based on measured signals and patient diary

Statement

Findings

Nearly all occurrences of respiratory events are in the supine position. The recording was technically successful and of high quality.

Conclusion

The recording findings are consistent with mild obstructive sleep apnea. Avoiding sleeping in a supine position (positional therapy), for example by placing a tennis ball in the back of the pajamas or wearing a sleep apnea belt, might be sufficient treatment. It is recommended to confirm the treatment outcome with a new recording.

23/02/2022

Dr. Sleep Apnea

This report is the consulting doctor's opinion based on the monitoring and the available background information. The final conclusions are always made by the patient's attending doctor who has the best

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Accept the report

Once you accept the report, the recording will be closed. Are you sure you want to accept the report?

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Patient ID: PID0a743_013381 Bittium HSAT

Close

Organization	Sleep test center	Created	13.12.2021 09:29	Last data transfer	13.12.2021 09:31
Status	● Ended Closed	Started	24.06.2021 00:00	Last update	29.12.2021 16:20
Assignment	Daniela DiagnosticUser	Ready for analysis	13.12.2021 09:44	Duration	1 / 1 day(s)

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Search: PID0a743 Create recording

10 All

Patient ID	Status	Last data transfer	Last update	Organization	Created	Type	Device	Assignment	Action
PID0a743	● Ended Closed	13.12.2021 09:31	29.12.2021 16:20	Sleep test center	13.12.2021 09:29	Bittium HSAT	0000552	Daniela DiagnosticUser	📄 ⬇️

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
Latest reports				
Recording	Organization	Added	Approval	
● PID0a74a	DemoClinic2	07.10.2020 17:49	07.10.2020 17:54	 

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Recordings

Search: PID0a74a Create recording - 10 - All -

Patient ID	Status	Last data transfer	Last update	Organization	Created	Type	Device	Assignment	Action
PID0a74a	● Closed	14.12.2021 15:47	17.12.2021 15:22	Sleep test center	14.12.2021 15:47	Bittium HSAT	0123456	None	 

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

Latest reports				
Recording	Organization	Updated	Approval	
● PID023d3	Sleep test center	16.12.2021 09:46 Not assigned	16.12.2021 09:48	 
● PID0a74a	Sleep test center	15.12.2021 15:16 Not assigned	Approve	 

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EN - Sales Technician Technician

Patient ID: PID0a74a Bittium HSAT Approve

Organization	Sleep test center	Created	14.12.2021 15:47	Last data transfer	14.12.2021 15:47
Status	● Closed	Started	24.02.2021 00:00	Last update	17.12.2021 15:22
Assignment	None	Ready for analysis	14.12.2021 15:49	Duration	1 / 1 day(s)

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Confirmation

Are you sure you want to approve the reports of Patient ID: 01-10-56? After approval, reports will be visible for the clinic.

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Patient ID: PID0a74a Bittium HSAT Approved

Organization	Sleep test center	Created	14.12.2021 15:47	Last data transfer	14.12.2021 15:47
Status	● Closed	Started	24.02.2021 00:00	Last update	31.12.2021 10:13
Assignment	None	Ready for analysis	14.12.2021 15:49	Duration	1 / 1 day(s)

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Figure 56 Enregistrement approuvé

Latest reports

Recording	Organization	Updated	Approval	
● PID023d3	Sleep test center	16.12.2021 09:46 Not assigned	16.12.2021 09:48	📄 ↓
● PID0a74a	Sleep test center	15.12.2021 15:16 Not assigned	16.12.2021 09:48	📄 ↓

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Latest reports

Recording	Organization	Updated	Approval	
● PID023d3	Sleep test center	16.12.2021 09:46 Not assigned	16.12.2021 09:48	 
● PID0a74a	Sleep test center	15.12.2021 15:16 Not assigned	16.12.2021 09:48	 

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The screenshot shows the Bittium MedicalSuite Center interface. At the top, there is a navigation bar with 'Back', 'Bittium', 'Dashboard', and 'Recordings'. The user is logged in as 'Sales Technician' in 'EN' language. The patient ID is 'PID0a870_013696' and the device is 'Bittium HSAT'. A table displays recording details:

Field	Value	Field	Value	Field	Value
Organization	Sleep test center	Created	21.12.2021 13:38	Last data transfer	21.12.2021 13:39
Status	Ended	Started	23.05.2021 22:53	Last update	23.12.2021 13:41
Assignment	None	Ready for analysis	21.12.2021 13:42	Duration	1 / 1 day(s)

Below the table, there are tabs for 'Recording data', 'Patient details', 'Anamneses', 'Patient diary', 'Recording details', and 'Attachments'. A 'Download' button is visible, along with a checkbox for 'Include all attachments'. At the bottom, there is a 'Recording data' table:

Day No.	Start date of file	Duration	Status	Quality	Notes
1/1	23.05.2021 19:53	08h 00m	Open	Poor	Add notes

A red box highlights the 'Close' and 'Reject' buttons in the top right corner of the recording details section.

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< Back **Bittium** Dashboard Recordings EN Sales Technician Technician

Patient ID: PID020c5_005914 Bittium HSAT

Organization	Evismo test center	Created	07.09.2021 15:04	Last data transfer	No date
Status	<input type="radio"/> Rejected	Started	No date	Last update	22.09.2021 14:37
Assignment	None	Ready for analysis	No date	Duration	

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Reject recording

Patient ID: **PID023d3**

The recorder 0000780 is currently in use in the selected recording. It will be released when the recording is rejected. Are you sure you want to reject the recording?

This action cannot be undone.

Why do you want to reject the recording?

- Corrupted recording data
- Not enough analyzable data due to noisy signal or unattached recorder
- Not enough sleep to write a statement
- Other

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Delete recording

Patient ID: **PID023d3**

The recorder 0000780 is currently in use in the selected recording. It will be released when the recording is deleted. Are you sure you want to delete the recording?

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The screenshot shows the Bittium user interface for a patient named PID0a74a. The top navigation bar includes 'Back', 'Bittium', 'Dashboard', 'Recordings', and user information 'EN - Sales Technician Technician'. The patient ID and name 'PID0a74a Bittium HSAT' are displayed, along with an 'Approve' button. Below this is a summary table of recording details:

Organization	Sleep test center	Created	16.12.2021 15:27	Last data transfer	16.12.2021 16:24
Status	Closed	Started	16.12.2021 15:36	Last update	21.12.2021 13:58
Assignment	None	Ready for analysis	17.12.2021 14:24	Duration	2 / 2 day(s)

Below the table are tabs for 'Recording data', 'Patient details', 'Anamneses', 'Patient diary', 'Recording details', and 'Attachments'. A 'Download' button is highlighted with a red box, along with an 'Include all attachments' checkbox. 'Cancel' and 'Save' buttons are also visible. The 'Recording data' section contains a table with the following data:

<input type="checkbox"/>	Day No.	Start date of file	Duration	Status	Quality	Notes	
<input type="checkbox"/>	1/2	16.12.2021 15:36	00h 29m	Open	Poor	Add notes	
<input checked="" type="checkbox"/>	2/2	17.12.2021 10:01	03h 00m	Open	Poor	Add notes	

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FABRICANT

Plateforme de services Bittium MedicalSuite

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Copyright 2024 Bittium Biosignals Ltd. Tous droits réservés. Les informations contenues dans ce document peuvent faire l'objet de modifications sans préavis. Bittium Biosignals Ltd conserve le droit de propriété et tous les autres droits sur le contenu de ce document. Toute reproduction du contenu de ce document est interdite sans l'accord écrit préalable de Bittium Biosignals Ltd.