# Bittium

---

## MedicalSuite™ – plattform HSAT-tjeneste

Klinisk veiledning

## Utgitt av

---

Bittium Biosignals Ltd.  
Pioneerinkatu 6  
70800 Kuopio  
Finland  
Telefon: +358 40 344 2000  
www.bittium.com

## Juridisk merknad

---

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
The screenshot shows the Bittium Recordings page. At the top, there is a navigation bar with the Bittium logo, 'Dashboard', and 'Recordings' tabs. On the right, it shows 'EN' and a user profile for 'Sales Technician'. Below the navigation bar, there is a search bar with a magnifying glass icon and the text 'Search', a 'Create recording' button, and two dropdown menus for '10' and 'All'. The main content is a table with the following columns: Patient ID, Status, Last data transfer, Last update, Organization, Created, Type, Device, Assignment, and Action. The table contains 11 rows of data, each representing a recording with various patient IDs, statuses (Active or Closed), dates, and device information.

Patient ID	Status	Last data transfer	Last update	Organization	Created	Type	Device	Assignment	Action
KK02	Active	28.10.2021 15:25	28.10.2021 15:46	Sleep sales center	28.10.2021 15:23	Bittium HSAT	0000568	Sales diag group 1	Download
Pat123456	Closed	19.10.2021 11:41	19.10.2021 12:12	Sleep sales center	19.10.2021 11:40	Bittium HSAT	0000673	Sales diag group 1	Download
PID0000b	Closed	18.10.2021 17:48	18.10.2021 18:16	Sleep sales center	18.10.2021 17:47	Bittium HSAT	0000673	Sales diag group 1	Download
PID0000a	Closed	18.10.2021 16:19	18.10.2021 16:15	Sleep sales center	18.10.2021 16:18	Bittium HSAT	0000536	Sales diag group 1	Download
PID00009	Active	16.10.2021 08:30	09.11.2021 10:16	Sleep sales center	16.10.2021 08:29	Bittium HSAT	0000536	Sales diag group 1	Download
PID00008	Active	16.10.2021 08:25	16.10.2021 09:04	Sleep sales center	16.10.2021 08:25	Bittium HSAT	0000568	Sales diag group 1	Download
Test008	Closed	15.10.2021 19:01	18.10.2021 15:26	Sleep sales center	15.10.2021 19:00	Bittium HSAT	0000673	Sales diag group 1	Download
TEST007	Closed	15.10.2021 18:46	15.10.2021 19:14	Sleep sales center	15.10.2021 18:43	Bittium HSAT	0000673	Sales diag group 1	Download
PID00007	Active	No date	14.10.2021 17:22	Sleep sales center	14.10.2021 17:18	Bittium HSAT	1234567	Sales diag group 1	Download
23_45_08	Closed	01.10.2021 16:17	14.10.2021 17:19	Sleep sales center	01.10.2021 16:16	Bittium HSAT	0123456	Sales diag group 1	Download

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The screenshot shows the Bittium Recordings page. At the top, there is a navigation bar with 'Bittium', 'Dashboard', and 'Recordings' tabs. The user is logged in as 'Sales Technician'. Below the navigation bar, there is a search bar and a 'Create recording' button. A table of recordings is displayed with columns: Patient ID, Status, Last data transfer, Last update, Organization, Created, Type, Device, and Assignment. A dropdown menu is open in the top right corner, showing options for the number of recordings to display: 10, 50, and 100. The number 10 is selected.

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The screenshot shows the Bittium Recordings page. At the top, there is a navigation bar with 'Bittium', 'Dashboard', and 'Recordings' tabs. The user is logged in as 'Sales Technician'. Below the navigation bar, there is a search bar and a 'Create recording' button. A table of recordings is displayed with columns: Patient ID, Status, Last data transfer, Last update, Organization, Created, Type, Device, and Assignment. A dropdown menu is open in the top right corner, showing options for the status of recordings: All, Active, and Closed. The 'All' option is selected.

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The screenshot shows the Bittium Recordings page. At the top, there is a navigation bar with 'Bittium', 'Dashboard', and 'Recordings' tabs. The user is logged in as 'Sales Technician'. Below the navigation bar, there is a search bar and a 'Create recording' button. A table of recordings is displayed with columns: Patient ID, Status, Last data transfer, Last update, Organization, Created, Type, Device, and Assignment. The 'Patient ID' column header has a dropdown arrow, and the 'Created' column header has a dropdown arrow. Both are highlighted with red boxes. The 'Patient ID' dropdown shows '22\_47\_08' and the 'Created' dropdown shows '11-10-2021' and '01-10-2021'.

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The screenshot shows the Bittium Recordings page. At the top, there is a navigation bar with 'Bittium', 'Dashboard', and 'Recordings' tabs. The user is logged in as 'Technician 1'. Below the navigation bar, there is a search bar and a 'Create recording' button. The main content is a table of recordings with columns: Patient ID, Status, Last data transfer, Last update, Organization, Created, Type, Device, Assignment, and Action. A dropdown menu is open for the 'Assignment' column of the first row, showing a search bar and a list of groups: 'Cardioscope Group', 'Diagnostic group', 'Diagnostic group (autoreports)', 'HRV group (autoreports)', and 'HSAT Analysis' (which is selected). Below the groups, there is a section for 'Users in HSAT Analysis' with 'Dr. Hanna HSAT' listed.

Patient ID	Status	Last data transfer	Last update	Organization	Created	Type	Device	Assignment	Action
testPDF	Closed	23.01.2024 15:14	25.01.2024 13:41	ServiceProvider	23.01.2024 15:14	Bittium	0000533	HSAT Analysis	[Download]
HSAT_ECG	Ended	16.01.2024 16:13	16.01.2024 16:17	DemoClinic1	16.01.2024 16:13	Bittium			[Download]
HSAT_2_RIP_Belts	Active	16.01.2024 16:11	25.01.2024 10:10	DemoClinic1	16.01.2024 16:11	Bittium			[Download]
PID00314	Ended	14.12.2023 15:06	19.12.2023 13:09	DemoClinic1	14.12.2023 15:05	Bittium			[Download]
PID002f3	Closed	22.11.2023 11:37	22.11.2023 11:41	DemoClinic1	22.11.2023 11:37	Bittium			[Download]

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The screenshot shows the Bittium Recordings page. At the top, there is a navigation bar with 'Bittium', 'Dashboard', and 'Recordings' tabs. The user is logged in as 'Sales Technician'. Below the navigation bar, there is a search bar and a 'Create recording' button. The main content is a table of recordings with columns: Patient ID, Status, Last data transfer, Last update, Organization, Created, Type, Device, Assignment, and Action. A dropdown menu is open for the 'Action' column of the first row, showing a search bar and a list of items: 'Recording data', 'Recording data 1 day(s)', 'Reports', '6169a6099ec5e8000acd5a02\_statement.pdf', and 'Other documents', 'recording-information.pdf'. A 'Download' button is visible at the top of the dropdown menu.

Patient ID	Status	Last data transfer	Last update	Organization	Created	Type	Device	Assignment	Action
Test008	Closed	15.10.2021 19:01	18.10.2021 15:26	Sleep sales center	15.10.2021 19:00	Bittium HSAT	0000533	Sales diag group 1	[Download]
TEST007	Closed	15.10.2021 18:46	15.10.2021 19:14	Sleep sales center	15.10.2021 18:43	Bittium HSAT	0000533		[Download]
PID0000b	Closed	18.10.2021 17:48	18.10.2021 18:16	Sleep sales center	18.10.2021 17:47	Bittium HSAT	0000533		[Download]
PID0000a	Closed	18.10.2021 16:19	18.10.2021 18:15	Sleep sales center	18.10.2021 16:18	Bittium HSAT	0000533		[Download]

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The screenshot displays the Bittium MedicalSuite Center interface. At the top, there is a navigation bar with a back arrow, the Bittium logo, and tabs for Dashboard and Recordings. The user is logged in as a Sales Technician. The main content area shows patient information for PID0a74a, including organization (Sleep sales center), status (Closed), and assignment (Sales diag group 1). A table of recording data is visible, with columns for Organization, Status, Assignment, Created, Started, Ready for analysis, Last data transfer, Last update, and Duration. Below the table are tabs for Recording data, Patient details, Anamneses, Patient diary, Recording details, and Attachments. The Patient details tab is active, showing a form with fields for Patient ID (Pat123456), Age (years) (55), Height (cm) (175), Weight (kg) (80), BMI (26.12), and Gender (Female). A note explains the Patient ID format: 'Patient ID must be at least 4 characters long and lower/uppercase letters, numbers and underscores are allowed. If Patient ID is left blank the system will generate it (for example PID0000A).' Buttons for Cancel and Save are located at the bottom right.

Organization	Sleep sales center	Created	19.10.2021 11:40	Last data transfer	19.10.2021 11:41
Status	● Closed	Started	18.10.2021 23:30	Last update	19.10.2021 12:12
Assignment	Sales diag group 1	Ready for analysis	No date	Duration	1 / 1 day(s)


Patient ID	Pat123456	Age (years)	55
<small>Patient ID must be at least 4 characters long and lower/uppercase letters, numbers and underscores are allowed. If Patient ID is left blank the system will generate it (for example PID0000A).</small>		Height (cm)	175
		Weight (kg)	80
		BMI	26.12
		Gender	Female

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

A screenshot of the Bittium MedicalSuite web application. The top navigation bar includes the Bittium logo, 'Dashboard', 'Recordings', and user information for 'Sales Technician'. The main content area is titled 'Create recording: Bittium HSAT' and has four tabs: 'Patient details\*', 'Anamneses', 'Recording details\*', and 'Attachments'. The 'Patient details' tab is active, showing a form with fields for Patient ID, Age (years)\*, Height (cm)\*, Weight (kg)\*, BMI\*, and Gender\*. A note explains that the Patient ID must be at least 4 characters long and contains instructions on how the system will generate it if left blank. Navigation buttons for 'Back', 'Next', 'Cancel', and 'Save' are visible.

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## Create recording: Bittium HSAT

Patient details Anamneses Recording details Attachments

<< Back Next >>

Cancel

Save

### Medication

Medication

Medication during the recording

No  Yes

### Medical history

Close family members diagnosed or treated for a sleep disorder

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## Chance of dozing off or falling asleep in the following situations

Sitting and reading	<input type="radio"/> Never	<input type="radio"/> Slight	<input type="radio"/> Moderate	<input type="radio"/> High
Watching TV	<input type="radio"/> Never	<input type="radio"/> Slight	<input type="radio"/> Moderate	<input type="radio"/> High
Sitting, inactive in a public place (e.g. a theatre or a meeting)	<input type="radio"/> Never	<input type="radio"/> Slight	<input type="radio"/> Moderate	<input type="radio"/> High
As a passenger in a car for an hour without a break	<input type="radio"/> Never	<input type="radio"/> Slight	<input type="radio"/> Moderate	<input type="radio"/> High
Lying down to rest in the afternoon when circumstances permit	<input type="radio"/> Never	<input type="radio"/> Slight	<input type="radio"/> Moderate	<input type="radio"/> High
Sitting and talking to someone	<input type="radio"/> Never	<input type="radio"/> Slight	<input type="radio"/> Moderate	<input type="radio"/> High
Sitting quietly after a lunch without alcohol	<input type="radio"/> Never	<input type="radio"/> Slight	<input type="radio"/> Moderate	<input type="radio"/> High
In a car, while stopped for a few minutes in the traffic	<input type="radio"/> Never	<input type="radio"/> Slight	<input type="radio"/> Moderate	<input type="radio"/> High

## Sleep history

Morning headache	<input type="radio"/> Never	<input type="radio"/> Rarely	<input type="radio"/> Sometimes	<input type="radio"/> Often	<input type="radio"/> Almost always
Dry mouth in the morning	<input type="radio"/> Never	<input type="radio"/> Rarely	<input type="radio"/> Sometimes	<input type="radio"/> Often	<input type="radio"/> Almost always
Daytime fatigue or snoring	<input type="radio"/> Never	<input type="radio"/> Rarely	<input type="radio"/> Sometimes	<input type="radio"/> Often	<input type="radio"/> Almost always
Daytime napping	<input type="radio"/> Never	<input type="radio"/> Rarely	<input type="radio"/> Sometimes	<input type="radio"/> Often	<input type="radio"/> Almost always
Heartburn	<input type="radio"/> Never	<input type="radio"/> Rarely	<input type="radio"/> Sometimes	<input type="radio"/> Often	<input type="radio"/> Almost always
Loss of concentration	<input type="radio"/> Never	<input type="radio"/> Rarely	<input type="radio"/> Sometimes	<input type="radio"/> Often	<input type="radio"/> Almost always
Memory loss	<input type="radio"/> Never	<input type="radio"/> Rarely	<input type="radio"/> Sometimes	<input type="radio"/> Often	<input type="radio"/> Almost always
Cramping at night	<input type="radio"/> Never	<input type="radio"/> Rarely	<input type="radio"/> Sometimes	<input type="radio"/> Often	<input type="radio"/> Almost always
Difficulties falling asleep	<input type="radio"/> Never	<input type="radio"/> Rarely	<input type="radio"/> Sometimes	<input type="radio"/> Often	<input type="radio"/> Almost always
Breathing breaks while sleeping	<input type="radio"/> Never	<input type="radio"/> Rarely	<input type="radio"/> Sometimes	<input type="radio"/> Often	<input type="radio"/> Almost always
Pain during night	<input type="radio"/> Never	<input type="radio"/> Rarely	<input type="radio"/> Sometimes	<input type="radio"/> Often	<input type="radio"/> Almost always
Snoring or snorting during sleep	<input type="radio"/> Never	<input type="radio"/> Rarely	<input type="radio"/> Sometimes	<input type="radio"/> Often	<input type="radio"/> Almost always
Waking up during the night or restless sleep	<input type="radio"/> Never	<input type="radio"/> Rarely	<input type="radio"/> Sometimes	<input type="radio"/> Often	<input type="radio"/> Almost always
A strong urge to move your feet or an uncomfortable sensation in your feet related to this urge during sleep	<input type="radio"/> Never	<input type="radio"/> Rarely	<input type="radio"/> Sometimes	<input type="radio"/> Often	<input type="radio"/> Almost always
Urinary urgency during sleep	<input type="radio"/> Never	<input type="radio"/> Rarely	<input type="radio"/> Sometimes	<input type="radio"/> Often	<input type="radio"/> Almost always
If yes, how many times per night?	<input type="radio"/> Once	<input type="radio"/> Twice	<input type="radio"/> More often		
Average amount of sleep	_____ Hours				

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**Create recording: Bittium HSAT**

Patient details\* Anamneses **Recording details** Attachments

<< Back Next >> Cancel Save

**Recording details**

Organization \_\_\_\_\_ Referring physician \_\_\_\_\_  
Requested report language English \_\_\_\_\_  
Reason for recording  
\_\_\_\_\_

**Device information**

Recorder S/N (21) 19 XXXXXXXX \_\_\_\_\_  
See the side of the device. **Please add numbers after "19"**.

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## Create recording: Bittium HSAT

Patient details Anamneses Recording details\* Attachments

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Patient ID: **PID0a74a** Bittium HSAT

Ready for analysis | ⋮

Organization	Sleep test center	Created	29.12.2021 12:06	Last data transfer	29.12.2021 12:20
Status	● Active -	Started	09.09.2021 23:30	Last update	29.12.2021 12:21
Assignment	None -	Ready for analysis	No date	Duration	2 / 2 day(s)

Recording data Patient details Anamneses **Patient diary** Recording details Attachments

Cancel Save

## Patient diary - 09.09.2021, 23:30 (06h 59m)

Please, ask the patient to answer the following questions about the recording.

When did you go to bed?	hh:mm	In your estimation, when did you fall asleep?	hh:mm
What time did you wake up in the morning?	hh:mm	In your estimation, how many hours did you sleep?	hh:mm
How did you sleep?		Did you get out of bed during the recording?	
How did you sleep compared to usual?			

Other comments

If the patient woke up during the recording, please write down the events.

+ Add new event

Event			
Time	00:51	Duration	hh:mm
Event	Recorder button pressed		

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
Status	Uploading	Created	01.04.2019 15:34	Last update	01.04.2019 15:34
Assignment	Mabel Hicks	Started	---	Duration	---

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
The screenshot displays the Bittium MedicalSuite interface. At the top, a metadata table shows details for 'Clinic A', including creation and update dates, and recording status. Below this, a navigation bar highlights 'Recording data'. A 'Download' section is visible with a 'Cancel' button and a 'Save' button. The main area is a table of recording data with columns for Day no., Start time, Duration, Status, Quality, and Notes. Two rows are shown, with the second row (Day 2/2) highlighted. A red box highlights the 'Preparing' status icon in the first column of the second row. To the right, a modal window titled 'UPLOADS' shows the patient ID and a list of upload progress items: 'Uploaded' (checked) and 'Preparing data - 2 / 2' (circular arrow icon).

Organization	Clinic A	Created	01.04.2019 15:34	Last data transfer	13.04.2019 12:00
Status	Preparing	Started	10.04.2019 12:00	Last update	13.04.2019 12:00
Assignment	None	Ended	No date	Duration	2 days

Day no.	Start time	Duration	Status	Quality	Notes
1 / 2	10.04.2019 22:00	08h 05m	Open	Fair	Add notes
2 / 2	11.04.2019 22:00	07h 50m	Open	Good	Add note

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The screenshot shows the Bittium patient interface. At the top, there is a navigation bar with 'Back', 'Bittium', 'Dashboard', and 'Recordings'. The user is identified as 'Sales Technician' in 'EN'. The patient ID is 'PID0a74a' and the device is 'Bittium HSAT'. A 'Ready for analysis' button is visible. Below this, a table displays patient details:

Organization	Sleep test center	Created	29.12.2021 11:22	Last data transfer	29.12.2021 11:32
Status	Active	Started	16.12.2021 02:44	Last update	29.12.2021 11:33
Assignment	None	Ready for analysis	No date	Duration	2 / 2 day(s)

Below the table are tabs for 'Recording data', 'Patient details', 'Anamneses', 'Patient diary', 'Recording details', and 'Attachments'. A 'Download' button is present with an option to 'Include all attachments'. A 'Recording data' section shows a table with one entry:

Day No.	Start date of file	Duration	Status	Quality	Notes
1/1	16.12.2021 00:44	08h 10m	Open	Good	Add notes

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Cancel Save

### Recording data

<input type="checkbox"/>	Day No.	Start date of file	Duration	Status	Quality	Notes
<input type="checkbox"/>	1/3	14.12.2021 23:42	04h 30m	Done	Fair	Add notes
<input type="checkbox"/>	2/3	14.12.2021 23:42	05h 00m	Open	Good	Add notes
<input type="checkbox"/>	3/3	14.12.2021 23:42	00h 15m	Unusable	Poor	Add notes

Marked as "Unusable"

- Corrupted recording data
- Not enough analyzable data due to noisy signal or unattached recorder
- Not enough sleep to write a statement
- Other (trs)

Figur 35 Velge status for opptaksdata

Unusable recording data

Patient ID: [patient\_id], 10/04/2019, 22:00 (08h 20m)

Why the recording data of this night is unusable?

Corrupted recording data

Not enough analyzable data due to noisy signal or unattached recorder

Not enough sleep to write the statement

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

The screenshot shows the patient details page for PID0a74a. At the top right, there is a 'Ready for analysis' button highlighted with a red box. Below the patient information, there is a table with columns for Organization, Status, Assignment, Created, Started, Ready for analysis, Last data transfer, Last update, and Duration. The 'Status' dropdown menu is highlighted with a red box, showing 'Active' and 'Ended' options. Below the table, there are tabs for Recording data, Patient details, Anamneses, Patient diary, Recording details, and Attachments. At the bottom, there is a 'Download' button and a 'Save' button.

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The screenshot shows the 'Recordings' page in the Bittium MedicalSuite Center. At the top, there is a search bar for PID and a 'Create recording' button. Below the search bar, there is a table with columns for Patient ID, Status, Last data transfer, Last update, Organization, Created, Type, Device, Assignment, and Action. The 'Status' dropdown menu for the first recording (PID023d3) is highlighted with a red box, showing 'Active' and 'Ended' options.

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## Ready for analysis

Patient ID **PID0a870\_013696**

Are you sure this recording is ready for analysis?

- Please make sure the patient diary is filled or it is attached to "Other documents" section
- You will no longer be able to upload new recording data to it.

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Diagnostic report needed

Recording	Assignment	Ready for analysis	Duration	
● PID0a74a	Test diag group 1	24.11.2021 09:10	1 day(s)	Analyze

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
The screenshot shows the Bittium interface for a patient record (Patient ID: PID00298). The top navigation bar includes 'Back', 'Bittium', 'Dashboard', 'Recordings', 'EN', a help icon, and the user profile 'Technician 1'. The patient record details include 'Organization: DemoClinic1', 'Created: 30.05.2023 13:56', 'Status: Active', 'Started: 23.11.2021 23:42', 'Ready for analysis: No date', and 'Duration: 1 / 1 day(s)'. A dropdown menu for 'Assignment' is open, showing a search bar and a list of groups: 'Cardiscope Group', 'Diagnostic group', 'Diagnostic group (autoreports)', 'HRV group (autoreports)', and 'HSAT Analysis' (which is selected). Below the groups, there is a section for 'Users in HSAT Analysis' with 'Dr. Hanna HSAT' listed. The interface also features a 'Ready for analysis' button, a 'Cancel' button, and a 'Save' button.

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## Diagnostic report needed

Recording	Assignment	Ready for analysis	Duration	
PatientZero	Test diag group 1	24.11.2021 09:10	1 day(s)	Analyze

Figur 44 Analyseknappen i instrumentpaneldialogboksen «Diagnoserapport nødvendig»

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Summary and statement | Analysis | Report

Accept report

Patient

Patient ID:	23-30-48	Height:	185cm	Weight:	85kg
BMI:	24.8	Age:	58 years	Gender:	Male

Summary

Apnea-hypopnea index: 8.7/h  
(Supine: 15.4/h, Non-supine: 2.1/h)

Oxygen Desaturation Index (ODI3): 8.7/h

AHI

Normal | Mild | Moderate | Severe

5 | 15 | 30

ODI3

Recording date	15/10/2021	Recording: TRT	08h 36m (23:30 - 08:07)
Evaluation period	08h 16m (23:43 - 08:00)		

Apnea-hypopnea index (AHI) displayed here is the total number of respiratory events (Apneas + Hypopneas) divided by Evaluation period. Evaluation period was determined based on measured signals and patient diary

Statement

Findings

Nearly all occurrences of respiratory events are in the supine position. The recording was technically successful and of high quality.

Conclusion

The recording findings are consistent with mild obstructive sleep apnea. Avoiding sleeping in a supine position (positional therapy), for example by placing a tennis ball in the back of the pajamas or wearing a sleep apnea belt, might be sufficient treatment. It is recommended to confirm the treatment outcome with a new recording.

23/02/2022

Dr. Sleep Apnea

This report is the consulting doctor's opinion based on the monitoring and the available background information. The final conclusions are always made by the patient's attending doctor who has the best

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## Accept the report

Once you accept the report, the recording will be closed. Are you sure you want to accept the report?

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Patient ID: PID0a743\_013381 Bittium HSAT

Close

Organization	Sleep test center	Created	13.12.2021 09:29	Last data transfer	13.12.2021 09:31
Status	Ended	Started	24.06.2021 00:00	Last update	29.12.2021 16:20
Assignment	Daniela DiagnosticUser	Ready for analysis	13.12.2021 09:44	Duration	1 / 1 day(s)

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Recordings

Search: PID0a743 Create recording

Patient ID	Status	Last data transfer	Last update	Organization	Created	Type	Device	Assignment	Action
PID0a743	Ended	13.12.2021 09:31	29.12.2021 16:20	Sleep test center	13.12.2021 09:29	Bittium HSAT	0000552	Daniela DiagnosticUser	Download

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
Latest reports				
Recording	Organization	Added	Approval	
● PID0a74a	DemoClinic2	07.10.2020 17:49	07.10.2020 17:54	 

Figur 51 Forhåndsvisningsikonet i dialogboksen «Siste rapporter»


Recordings									
Q PID0a74a		Create recording -		10 -		All -			
Patient ID	Status	Last data transfer	Last update	Organization	Created	Type	Device	Assignment	Action
PID0a74a	● Closed	14.12.2021 15:47	17.12.2021 15:22	Sleep test center	14.12.2021 15:47	Bittium HSAT	0123456	None	 

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Latest reports				
Recording	Organization	Updated	Approval	
● PID023d3	Sleep test center	16.12.2021 09:46 Not assigned	16.12.2021 09:48	 
● PID0a74a	Sleep test center	15.12.2021 15:16 Not assigned	<b>Approve</b>	 

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The screenshot shows the Bittium interface with a patient record for PID0a74a. The 'Approve' button is highlighted with a red box. Below the patient details is a navigation bar with tabs for Recording data, Patient details, Anamneses, Patient diary, Recording details, and Attachments. At the bottom, there are download options and 'Cancel' and 'Save' buttons.

Organization	Sleep test center	Created	14.12.2021 15:47	Last data transfer	14.12.2021 15:47
Status	<span style="color: red;">●</span> Closed	Started	24.02.2021 00:00	Last update	17.12.2021 15:22
Assignment	None	Ready for analysis	14.12.2021 15:49	Duration	1 / 1 day(s)

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The confirmation dialog box has a blue header with the word 'Confirmation'. The main text asks: 'Are you sure you want to approve the reports of Patient ID: 01-10-56? After approval, reports will be visible for the clinic.' At the bottom, there are 'Cancel' and 'Approve' buttons.

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< Back **Bittium** Dashboard Recordings EN Sales Technician Technician

**Patient ID: PID0a74a** Bittium HSAT Approved

Organization	Sleep test center	Created	14.12.2021 15:47	Last data transfer	14.12.2021 15:47
Status	<span style="color: red;">●</span> Closed	Started	24.02.2021 00:00	Last update	31.12.2021 10:13
Assignment	None	Ready for analysis	14.12.2021 15:49	Duration	1 / 1 day(s)

Recording data Patient details Anamneses Patient diary Recording details Attachments

Download Cancel Save

Include all attachments

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Latest reports

Recording	Organization	Updated	Approval
<span style="color: red;">●</span> PID023d3	Sleep test center	16.12.2021 09:46 Not assigned	16.12.2021 09:48
<span style="color: red;">●</span> PID0a74a	Sleep test center	15.12.2021 15:16 Not assigned	16.12.2021 09:48

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
## Latest reports

Recording	Organization	Updated	Approval	
● PID023d3	Sleep test center	16.12.2021 09:46 Not assigned	16.12.2021 09:48	 
● PID0a74a	Sleep test center	15.12.2021 15:16 Not assigned	16.12.2021 09:48	 

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3 file(s)

+ Add file or drag here

<input type="checkbox"/>	Filename	Added			
<input type="checkbox"/>	recording-information.pdf				
<input type="checkbox"/>	Document_1.txt	19.12.2021 14:53			

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Navigation: < Back | Bittium | Dashboard | Recordings | EN | Sales Technician | Technician

Patient ID: PID0a870\_013696 Bittium HSAT

Organization	Sleep test center	Created	21.12.2021 13:38	Last data transfer	21.12.2021 13:39
Status	● Ended	Started	23.05.2021 22:53	Last update	23.12.2021 13:41
Assignment	None	Ready for analysis	21.12.2021 13:42	Duration	1 / 1 day(s)

Recording data | Patient details | Anamneses | Patient diary | Recording details | Attachments

Download  Include all attachments [Cancel] [Save]

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Navigation: < Back | Bittium | Dashboard | Recordings | EN | Sales Technician | Technician

Patient ID: PID020c5\_005914 Bittium HSAT

Organization	Evismo test center	Created	07.09.2021 15:04	Last data transfer	No date
Status	○ Rejected	Started	No date	Last update	22.09.2021 14:37
Assignment	None	Ready for analysis	No date	Duration	

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## Reject recording

Patient ID: **PID023d3**

The recorder 0000780 is currently in use in the selected recording. It will be released when the recording is rejected. Are you sure you want to reject the recording?

**This action cannot be undone.**

Why do you want to reject the recording?

- Corrupted recording data
- Not enough analyzable data due to noisy signal or unattached recorder
- Not enough sleep to write a statement
- Other

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## Delete recording

Patient ID: **PID023d3**

The recorder 0000780 is currently in use in the selected recording. It will be released when the recording is deleted. Are you sure you want to delete the recording?

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The screenshot shows the Bittium Recordings page. The top navigation bar includes 'Bittium', 'Dashboard', and 'Recordings'. The user is logged in as 'Sales Technician'. A search bar is present with the text 'Data'. A table of recordings is displayed with columns: Patient ID, Last data transfer, Last update, Organization, Created, Type, Device, Assignment, and Action. A tooltip is shown over a recording with the status 'Rejected', listing reasons: 'Not enough analyzable data due to noisy signal or unattached recorder' and 'Not enough sleep during the recording to write the statement'.

Patient ID	Last data transfer	Last update	Organization	Created	Type	Device	Assignment	Action
DataTransferTest	20.12.2021 14:58	21.12.2021 15:11	Sleep test center	20.12.2021 14:06	Bittium HSAT Plus	0000516	None	↓

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The screenshot shows the Bittium single recording page. The top navigation bar includes 'Bittium', 'Dashboard', and 'Recordings'. The user is logged in as 'Sales Technician'. The page displays details for a recording, including Patient ID, Organization, Status, Created, Started, Last data transfer, and Last update. A tooltip is shown over the 'Rejected' status, listing reasons: 'Not enough analyzable data due to noisy signal or unattached recorder' and 'Not enough sleep during the recording to write the statement'.

Organization	Created	Last data transfer
Bittium HSAT	20.12.2021 14:06	20.12.2021 14:58

Status	Started	Last update
Rejected	20.12.2021 14:09	21.12.2021 15:11

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Navigation: < Back | Bittium | Dashboard | Recordings | EN - | Sales Technician | Technician

Patient ID: PID0a74a | Bittium HSAT | Approve

Organization	Sleep test center	Created	16.12.2021 15:27	Last data transfer	16.12.2021 16:24
Status	<span style="color: red;">●</span> Closed	Started	16.12.2021 15:36	Last update	21.12.2021 13:58
Assignment	None	Ready for analysis	17.12.2021 14:24	Duration	2 / 2 day(s)

Recording data | Patient details | Anamneses | Patient diary | Recording details | Attachments

Download  
 Include all attachments

Cancel Save

### Recording data

<input type="checkbox"/>	Day No.	Start date of file	Duration	Status	Quality	Notes	
<input type="checkbox"/>	1/2	16.12.2021 15:36	00h 29m	Open	Poor	Add notes	
<input checked="" type="checkbox"/>	2/2	17.12.2021 10:01	03h 00m	Open	Poor	Add notes	

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