

Productive Business Management With Sophisticated Tools

The user and device management, as well as reporting are centrally managed with the MedicalSuite™ Back-Office tool and all information is securely accessible, regardless of location. In addition to creating and managing accounts, the Business User account offers accurate insights on the hospital's usage of the service. Customer statistics are immediately available for forecasting, reporting and device stock management.



More Information of Bittium MedicalSuite™



Scan to visit Bittium MedicalSuite™ website

Bittium

Connectivity
to be trusted.

Bittium is a trusted Nordic company with over 30 years of expertise in advanced biosignal processing. Bittium provides medical technology in the measurement of biosignals and monitoring for cardiology, neurophysiology, neuroscience, rehabilitation, occupational health and sports medicine. The products meet European Union medical CE requirements and the company's quality system meets ISO 9001 and ISO 13485 directive MDD 93/42/EEC requirements.

Cardiology

Bittium Faros™ 180
Bittium Faros™ 360
Bittium Cardiac Navigator™ Software

Neurophysiology and Neuroscience

Bittium BrainStatus™
Bittium NeurOne™

Sleep Apnea

Bittium Respiro™***
Bittium Respiro Analyst™ Software**

MedicalSuite™ Service Products

Bittium Holter™**
Bittium HolterPlus™**
Bittium HSAT™**
Bittium HSAT Plus™**

**A system that utilizes CE certified Bittium Faros™ and Bittium Cardiac Navigator™ products.*

***Bittium Respiro™ technology is not yet available for commercial use.*

CE certification is expected in H1/2022. and Bittium Cardiac Navigator™ products.

FOR MORE INFORMATION, PLEASE CONTACT:

medical@bittium.com
www.bittium.com

Bittium MedicalSuite™

Remote Monitoring Service Platform

For a Higher Level of Collaborative Care



Copyright © 2021-11-04 Bittium. All rights reserved.

Bittium

Efficient Collaboration Through the Entire Treatment Workflow

Bittium MedicalSuite™ is an advanced web-based remote monitoring service platform that facilitates secure information sharing between service providers, clinics and cardiologists to enable better collaborative care. The Service is a combination of intelligent analysis, wireless data transfer, reporting and management modules which can be optimized according to your needs. As a result, you can benefit from effortless workflows and enhanced efficiency.

Remote Monitoring Made Easy

Our outstanding capability to provide solutions for remote monitoring of patients outside the hospital creates new opportunities for minimizing healthcare costs. With Bittium MedicalSuite™ service platform the analysis and diagnostics services can be offered remotely from one digital channel to several locations. MedicalSuite™ service platform supports cardiac monitoring, Home Sleep Apnea Testing (HSAT) and HRV analysis.

One Platform - Multiple Applications

Cardiac Monitoring

Cardiac monitoring solution enables possibilities for traditional Holter recordings and long-term cardiac monitoring up to 30 days. Waterproof Bittium Faros™ ECG device, Bittium OmegaSnap™ ECG Electrodes and MedicalSuite platform flexibility make it possible to select the right solution based on your needs.

Bittium MedicalSuite™


Enabling better collaborative care with high-quality remote diagnostics

Diagnostics



Viewing your own work queue
Daily checks for meaningful events
Analysis and report

Service Center



Easy verification for received recording quality
Assign recordings for analysis to the specialists
Publishing diagnostics reports to clinics

Clinics



Setting up a recording
Uploading data for analysis
Receiving and downloading of reports

HRV Analysis

The analysis bases on the heart rate variability (HRV) which reflects the functioning of the autonomic nervous system. Bittium Faros™ ECG recording device combined with Bittium HRV Analyst™ software provides reliable data for proper analysis of overall health and well-being as part of Bittium MedicalSuite™ platform offering.

Home Sleep Apnea Testing

Bittium MedicalSuite™ provides the Home Sleep Apnea Testing™ (HSAT) services. HSAT solution is intended to record and analyze typical sleep related breathing disorders, such as obstructive sleep apnea, central sleep apnea, mixed sleep apnea, hypopnea and Cheyne-Stokes breathing. HSAT solution is combination of Bittium Respiro™ screening device, analysis software and MedicalSuite™ service platform.

